



CAPSLink 2.0.7 Release Notes

Cyclic Rate Orders

- Cyclic rate orders will now display the correct Order Info on the CAPSLink label.

Edit / Freeze issue

- When editing an existing previous day's order the system should no longer fail to re-open the Order Screen, which produced a "locked state" only recoverable with a refresh of the browser.

Order Screen Navigation Bar

- The Number of Units field on the Order Screen was moved to the right of 'print' and is more noticeable. The Navigation bar containing the Number of Units field is now "sticky", so when a user scrolls through the order it is always visible.

Prescriber List Sorting

- The Prescriber list now sorts in ascending order regardless of the letter case (upper or lower).

MRN Duplicates

- When entering a duplicate MRN for a new patient, additional instructions will appear on the resulting prompt to the user.

Ingredient Quantity Text Highlight

- Clicking in the ingredient quantity field on the Order Screen will now highlight the text, rather than placing the cursor to the right of the text.

Validated Order / View Errors Window

- The view errors window for an already validated order will no longer display override messages to the user.

Order List Quantities

- Quantities have been added to the radio buttons on the Order List for General, Time, and Status.